iVascular® sergeant

EN	INSTRUCTIONS FOR USE. PERIPHERAL SUPPORT CATHETER	03
ES	INSTRUCCIONES DE USO. CATÉTER DE SOPORTE PERIFÉRICO	03
PT	INSTRUÇÕES DE UTILIZAÇÃO. CATETER DE SUPORTE PERIFÉRICO	04
IT	ISTRUZIONI PER L'USO. CATETERE DI SUPPORTO PERIFERICO	05
FR	MODE D'EMPLOI. CATHÉTER DE SUPPORT PÉRIPHÉRIQUE	06
DE	GEBRAUCHSANWEISUNG. PERIPHERER SUPPORTKATHETER	06
NL	GEBRUIKSAANWIJZING. PERIFERE ONDERSTEUNINGSKATHETER	07

ENGLISH

1. Description of the product

The peripheral support catheter is a catheter with a single lumen from the proximal end to the distal end, and it is used both for the passage of a guidewire and for the delivery of contrast media or saline solution.

The proximal end of the catheter houses a female luer-lock connector, which can be coupled to male luer-lock accessories such as infusion syringes or haemostasis valves.

The distal part ends in a 5 mm long, rounded and atraumatic radiopaque tip to prevent damage to arteries as it advances. It may be straight or slightly angled.

From the tip, spaced 4 cm apart, it has three 4 mm radiopaque markers, which help when estimating the length and geometry of the area of the artery to be treated.

The distal body of the catheter has a hydrophilic coating, applied by immersion and subsequent UV curing. It is durable and helps with catheter navigation through tortuous arteries.

The useful length of the catheter may be 65 cm, 90 cm, 130 cm or 150 cm depending on the reference.

The product can have a straight or curved finish. The product may be equipped with guides of 0.014 inches, 0.018 inches or 0.035 inches.

Contents

- One peripheral support catheter inserted in a dispenser to avoid damage to the catheter, inside a sterile bag.
- A leaflet with the instructions for use.

2. Intended use / indication

The catheter is intended to facilitate, guide and support a guidewire during access into the peripheral vasculature, allow for guidewire exchanges, and/or provide a conduit for the delivery of saline solutions and/or diagnostic contrast agents.

The catheter is indicated for patients with peripheral artery disease in small vessel or challenging anatomy during diagnostic and interventional procedures.

3. Contraindications

Patients for whom antiplatelet and/or anticoagulant therapy is contraindicated.

4. Warnings

- The device must be used by experienced physicians familiar with the technique of percutaneous transluminal angioplasty (PTA).
- The percutaneous intervention must be carried out in hospitals equipped with emergency resources to develop open surgery or in hospitals with easy access to other hospitals that can perform this surgery.
- This is a single-use product. It must not be resterilised or reused once the procedure has finished. Reusing the product in another patient can cause cross-contamination, infections or the transmission of infectious diseases from one patient to another. Reusing the product can alter it and limit its efficacy.
- Reusing, reprocessing or resterilising can compromise the structural integrity of the device and/or cause it to fail, which at the same time can result in injury to the patient, disease or death.
- The product is supplied sterile. Check the expiry date and do not use products that have passed this date.
- Inspect the packaging before opening. If you find any defect or the packaging is damaged, discard the product and return it to the manufacturer.
- Use aseptic techniques when extracting the product from the primary packaging.
- Before using the catheter, inspect it to verify it is not curved or bent and dispose of any catheters that are not in good condition.
- Do not dry with gauze.
- Do not expose the product to organic solvents, or to contrast media that is not indicated for intravascular use.
- Proceed with extreme care to avoid damaging the catheter as it advances.
- Once inside the vascular system, handle the catheter under high-quality fluoroscopic vision. Do not allow the product to advance without the guidewire inside it. If you note any resistance, determine the cause

before continuing with the procedure or remove the catheter.

- When changing a guide with the product inside the vessel, take extreme care. If, when you insert the new guide, you notice resistance or that the catheter is bent, extract the product and the guide together and check that they are in good condition.
- Avoid introducing air or gaseous media through the inside of the catheter.
- Do not exceed the maximum injection pressure recommended. Follow the flow tables if you are performing an infusion through the catheter.
- This catheter is not indicated for neurovasculature or for coronary arteries.

5. Precautions

- Administer the appropriate medical treatment to the patient: anticoagulants, vasodilators, etc., according to the intravascular catheter insertion protocol.
- Before use, check the size of the catheter and its suitability for the intended use.
- Proceed with extreme care to avoid damaging the catheter as it advances. Take great care not damage the catheter in advance. If you notice any resistance as it advances, stop and determine the cause before continuing.
- If you notice any resistance when withdrawing the catheter, remove the catheter, the guidewire and the introducer all together.
- Store in a cool and dry place, away from direct sunlight.
- After use, this product can be a biological hazard. Handle and dispose of it in accordance with accepted medical practices, and the relevant local, state and federal laws and standards.

6. Possible adverse effects/complications

The possible adverse effects and/or complications that may become evident before, during or after the procedure include:

- Death
- Stroke/embolism/thrombosis
- Haemodynamic deterioration
- Arrhythmias
- Infections
- Total occlusion of the artery
- Perforation or dissection of the treated area
- Aneurysm or pseudoaneurysm
- Spasms
- Arteriovenous fistula
- Local haemorrhage with haematoma in the access area
- Hypo/hypertension
- Allergy to the materials

7. Instructions for use

7.1. Necessary equipment

- Heparinised normal saline solution.
- Contrast medium (use contrast media indicated for intravascular use).
- Guidewire of the size indicated on the label.
- Introducer with haemostatic valve of the size indicated on the label. Do not use smaller sizes that could damage the catheter.
- Various standard syringes of 10-20 cc with saline solution for cleaning the system.

7.2. Preparation of the catheter

- Remove the catheter from the protective dispenser. Inspect the catheter, if you notice any anomaly, not use it. Check that the size is suitable.
- Connect a 10 ml syringe containing sterile saline solution to the proximal luer-lock connector, and irrigate the lumen until the liquid comes out of the tip.
- To activate the hydrophilic coating and improve the catheter's navigability, dip the catheter in saline solution before use.

7.3. Insertion technique/treatment

- Fit the guide inside the patient according to PTA techniques, using fluoroscopy procedures to determine the position at all times.

- Insert the proximal end of the guide into the distal tip of the support catheter.
- Carefully advance the support catheter on the guidewire, through the introducer, aided by fluoroscopy. Ensure that the guide exits through the proximal luer connector of the support catheter.
- Advance the support catheter to the desired point with the help of fluoroscopic vision. Advance the guide with the catheter providing support and helping it pass the lesion.
- Once the goal has been achieved:
 - You can change the guide if necessary, extracting the first and placing the second.
 - You can inject a contrast liquid to see the vascular bed. To do so, extract the guidewire and connect the syringe to the proximal luer connector of the catheter.
- Once it reaches the lesion and when the guide has passed it, you can remove the support catheter.

8. Warranty

The product and all of its components have been designed, manufactured, tested, and packaged with maximum precaution. LVD Biotech guarantees the product until its expiry date provided the packaging is not broken, manipulated or damaged.

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